

Anti-Hnmt Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P50135

Item Number: CM0017698



Introduction

High-quality rabbit polyclonal anti-HNMT antibody validated for WB, IP, ELISA, targeting human and mouse histamine N-methyltransferase (33kDa). Ideal for histamine metabolism and airway response research.

[Learn More](#)

Parameter	Value
Product Name	HNMT Rabbit pAb
Remarks/Alias	HMT; MRT51; HNMT-S1; HNMT-S2; HNMT
Species	Human
GeneID	3176
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-292 of human HNMT (NP_008826.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P50135
Protein Weight	33kDa
Shipping	Ice bag