


Anti-Hnrnp A1 Polyclonal Antibody For Wb, If/Icc, Elisa - P09651

Item Number: CM0021497



Introduction

Rabbit polyclonal antibody against human hnRNP A1, validated by KO, for Western blot, immunofluorescence/ICC, and ELISA. Cross-reacts with mouse and rat. Research use for mRNA processing, splicing, and transport studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] hnRNP A1 Rabbit pAb
Remarks/Alias	UP 1; ALS19; ALS20; HNRPA1; IBMPFD3; HNRPA1L3; hnRNP A1; hnRNP-A1; A1
Species	Human
Gene ID (Human)	3178
Gene ID	3178
Immunogen	Recombinant protein containing aa 1-140 of human hnRNP A1 (NP_002127.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09651
Protein Weight	39kDa
Shipping	Ice bag