


Anti-Hnrnp A1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P09651

Item Number: CM0025611



Introduction

Rabbit polyclonal antibody targeting human hnRNP A1 (P09651). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with Human, Mouse, Rat. Ideal for RNA splicing, transport, and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	hnRNP A1 Rabbit pAb
Remarks/Alias	UP 1; ALS19; ALS20; HNRPA1; IBMPFD3; HNRPA1L3; hnRNP A1; hnRNP-A1
Species	Human
Gene ID (Human)	3178
Gene ID	3178
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-372 of human hnRNP A1 (NP_112420.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09651
Protein Weight	39kDa
Shipping	Ice bag