


Anti-Hnrnp A1 Rabbit Monoclonal Antibody For Wb - P09651

Item Number: CM0003158



Introduction

High-quality rabbit monoclonal antibody against human hnRNP A1 (P09651). Validated for WB, IHC-P, IF/ICC, IP, and ELISA. Cross-reacts with mouse and rat. Ideal for studying mRNA processing and neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	hnRNP A1 Rabbit mAb
Remarks/Alias	UP 1; ALS19; ALS20; HNRPA1; IBMPFD3; HNRPA1L3; hnRNP A1; hnRNP-A1
Species	Human
GeneID (Human)	3178
GeneID	3178
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 250-350 of human hnRNP A1 (P09651)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09651
Protein Weight	39kDa
Shipping	Ice bag