

Anti-Hnrnp A2B1 Rabbit Monoclonal Antibody [Kd Validated] For Wb, Ihc-P, If/Icc, Elisa - P22626

Item Number: CM0008190



Introduction

Anti-hnRNP A2B1 rabbit monoclonal antibody, KD validated, for Western blot, IHC, IF, and ELISA. Detects human, mouse, and rat hnRNP A2B1. Peptide immunogen; 37 kDa. Ideal for RNA processing and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] hnRNP A2B1 Rabbit mAb
Remarks/Alias	RNPA2; HNRPA2; HNRPB1; SNRPB1; HNRNPA2; HNRNPB1; IBMPFD2; HNRPA2B1
Species	Human
GeneID (Human)	3181
GeneID	3181
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 151-250 of human HNRNPA2B1 (NP_112533.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22626
Protein Weight	37kDa
Shipping	Ice bag