

Anti-Hnrnp C Rabbit Monoclonal Antibody For Wb, If, Elisa - P07910

Item Number: CM0010863



Introduction

Anti-hnRNP C rabbit monoclonal antibody, reactive with human and mouse, validated for WB, IF/ICC, IF-P, and ELISA. Targets heterogeneous nuclear ribonucleoproteins C1/C2, involved in pre-mRNA binding and splicing.

[Learn More](#)

Parameter	Value
Product Name	hnRNP C Rabbit mAb
Remarks/Alias	C1; C2; HNRNP; HNRPC; SNRPC; hnRNP C
Species	Human
GeneID (Human)	3183
GeneID	3183
Immunogen	Recombinant protein containing a sequence corresponding to amino acids 80-170 of human hnRNP C (P07910)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P07910
Protein Weight	34kDa
Shipping	Ice bag