

Anti-Hnrnp C1/C2 Rabbit Polyclonal Antibody - P07910

Item Number: CM0019847



Introduction

Rabbit polyclonal antibody against human hnRNP C (C1/C2). Validated for WB, IHC-P, IF/ICC, and ELISA; reacts with human, mouse, rat. Ideal for studies on pre-mRNA binding, mRNA stability, and splicing.

[Learn More](#)

Parameter	Value
Product Name	hnRNP C Rabbit pAb
Remarks/Alias	C1; C2; HNRNP; HNRPC; SNRPC; hnRNP C
Species	Human
GeneID (human)	3183
GeneID	3183
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 86-190 of human hnRNP C (NP_112604.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07910
Protein Weight	34kDa
Shipping	Ice bag