

## Anti-Hnrnp E1/Pcbp1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15365

Item Number: CM0017171



### Introduction

Human hnRNP E1/PCBP1 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat. Targets a single-stranded nucleic acid binding protein involved in RNA processing and viral replication. Ideal for IVD assay development.

[Learn More](#)

| Parameter        | Value                                                                                                                                        |
|------------------|----------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | hnRNP E1/PCBP1 Rabbit pAb                                                                                                                    |
| Remarks/Alias    | HNRPX; HNRPE1; hnRNP-X; HEL-S-85; hnRNP-E1; hnRNP E1/PCBP1                                                                                   |
| Species          | Human                                                                                                                                        |
| GeneID (human)   | 5093                                                                                                                                         |
| GeneID           | 5093                                                                                                                                         |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human hnRNP E1/hnRNP E1/PCBP1 (NP_006187.2). |
| Source           | Rabbit                                                                                                                                       |
| Category         | Polyclonal Antibodies                                                                                                                        |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                                                     |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                            |
| SWISS            | Q15365                                                                                                                                       |
| Protein Weight   | 37kDa                                                                                                                                        |
| Shipping         | Ice bag                                                                                                                                      |