

Anti-Hnrnp E1/Pcbp1 Rabbit Monoclonal Antibody For Wb, Elisa - Q15365

Item Number: CM0003667



Introduction

Recombinant rabbit monoclonal antibody targeting human PCBP1 (hnRNP E1). Validated for WB and ELISA; cross-reacts with mouse and rat. Poly(rC)-binding protein involved in RNA processing and erythropoiesis.

[Learn More](#)

Parameter	Description
Product Name	hnRNP E1/PCBP1 Rabbit mAb
Remarks/Alias	HNRPX; HNRPE1; hnRNP-X; HEL-S-85; hnRNP-E1; hnRNP E1/PCBP1
Species	Human
GeneID (Human)	5093
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PCBP1 (Q15365)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15365
Protein Weight	37kDa
Shipping	Ice bag