

Anti-Hnrnp E1/Pcbp1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q15365

Item Number: CM0004025



Introduction

High-performance rabbit monoclonal antibody to human hnRNP E1/PCBP1 (SWISS: Q15365). Validated for WB, IF/ICC, ELISA with cross-reactivity to human and mouse. Ideal for RNA-binding protein and erythropoiesis research.

[Learn More](#)

Parameter	Value
Product Name	hnRNP E1/PCBP1 Rabbit mAb
Remarks/Alias	HNRPX; HNRPE1; hnRNP-X; HEL-S-85; hnRNP-E1; hnRNP E1/PCBP1
Species	Human
GeneID	5093
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human hnRNP E1/PCBP1 (NP_006187.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15365
Protein Weight	37kDa
Shipping	Ice bag