


Anti-Hnrnp E2/Pcbp2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15366

Item Number: CM0024704



Introduction

Rabbit polyclonal antibody targeting human hnRNP E2/PCBP2 (Q15366), validated in WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying RNA processing and antiviral innate immunity.

[Learn More](#)

Parameter	Value
Product Name	hnRNP E2/PCBP2 Rabbit pAb
Remarks/Alias	HNRPE2; HNRNPE2; hnRNP-E2; hnRNP E2/PCBP2
Species	Human
GeneID (Human)	5094
GeneID	5094
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 150-250 of human hnRNP E2/PCBP2 (NP_001122383.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15366
Protein Weight	39kDa
Shipping	Ice bag