


Anti-Hnrnp E2/Pcbp2 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q15366

Item Number: CM0005415



Introduction

Rabbit monoclonal antibody targeting hnRNP E2/PCBP2 (SWISS: Q15366). Suitable for WB, IP, ELISA with human, mouse, rat reactivity. For IVD research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	hnRNP E2/PCBP2 Rabbit mAb
Remarks/Alias	HNRPE2; HNRNPE2; hnRNP-E2; hnRNP E2/PCBP2
Species	Human
Gene ID (Human)	5094
Gene ID	5094
Immunogen	Synthetic peptide - A synthetic peptide corresponding to a sequence within amino acids 150-250 of human hnRNP E2/PCBP2 (NP_001122383.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15366
Protein Weight	39kDa
Shipping	Ice bag