

Anti-Hnrrnp K Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P61978

Item Number: CM0022509



Introduction

Rabbit polyclonal antibody targeting human hnRNP K (P61978), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Recombinant immunogen (aa 1-280). Ideal for studies on pre-mRNA metabolism and DNA damage response.

[Learn More](#)

Parameter	Value
Product Name	hnRNP K Rabbit pAb
Remarks/Alias	AUKS; CSBP; TUNP; HNRPK; hnRNP K
Species	Human
GeneID (human)	3190
GeneID	3190
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human hnRNP K (NP_112552.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61978
Protein Weight	51kDa
Shipping	Ice bag