

Anti-Hnrnp K Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P61978

Item Number: CM0002689



Introduction

Rabbit monoclonal antibody targeting human hnRNP K; validated for Western blot, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for pre-mRNA processing, p53 signaling, and neuronal differentiation research.

[Learn More](#)

Parameter	Value
Product Name	hnRNP K Rabbit mAb
Remarks/Alias	AUKS; CSBP; TUNP; HNRPK; hnRNP K
Species	Human
GeneID (Human)	3190
GeneID	3190
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human hnRNP K (P61978).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61978
Protein Weight	49kDa/51kDa
Shipping	Ice bag