


Anti-Hnrnp Q/Syncrip Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O60506

Item Number: CM0012011



Introduction

Rabbit monoclonal antibody targeting human hnRNP Q/SYCRIP (SWISS: O60506). Validated for WB, IHC-P, IF/ICC, and ELISA; cross-reacts with mouse and rat. Suitable for RNA metabolism and translational control research.

[Learn More](#)

Parameter	Value
Product Name	hnRNP Q/SYCRIP Rabbit mAb
Remarks/Alias	PP68; NSAP1; GRYRBP; HNRNPQ; HNRPQ1; GRY-RBP; hnRNP-Q; hnRNP Q/SYCRIP
Species	Human
GeneID (Human)	10492
GeneID	10492
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 524-623 of human hnRNP Q/SYCRIP (O60506)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60506
Protein Weight	70kDa
Shipping	Ice bag