


Anti-Hnrnp Q/Syncrip Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O60506

Item Number: CM0030908



Introduction

Rabbit polyclonal antibody targeting hnRNP Q/SYCRIP, validated for WB, IHC-P, ELISA applications. Detects endogenous hnRNP Q in human, mouse, and rat samples. Ideal for mRNA processing and translation regulation studies.

[Learn More](#)

Parameter	Value
Product Name	hnRNP Q/SYCRIP Rabbit pAb
Remarks/Alias	PP68; NSAP1; GRYRBP; HNRNPQ; HNRPQ1; GRY-RBP; hnRNP-Q; hnRNP Q/SYCRIP
Species	Human
GeneID (Human)	10492
GeneID	10492
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human hnRNP Q /hnRNP Q/SYCRIP (NP_001153146.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60506
Protein Weight	70kDa
Shipping	Ice bag