

Anti-Hnrnp U Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q00839

Item Number: CM0028031

Introduction

Rabbit polyclonal antibody targeting hnRNP U (SAF-A), validated for WB, IHC-P, IF/ICC, IP, and ELISA in human, mouse, and rat. Suitable for chromatin architecture, transcription, and RNA biology studies.

[Learn More](#)

Parameter	Value
Product Name	hnRNP U Rabbit pAb
Remarks/Alias	SAFA; DEE54; HNRPU; SAF-A; U21.1; pp120; EIEE54; GRIP120; hnRNP U; HNRNPU-AS1
Species	Human
Gene ID (Human)	3192
Gene ID	3192
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 608-825 of human hnRNP U (NP_114032.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00839
Protein Weight	91kDa
Shipping	Ice bag