

Anti-Hnrnpa0 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q13151

Item Number: CM0028975

Introduction

Rabbit polyclonal antibody against human HNRNPA0 (HNRPA0). Validated for WB, IP, ELISA applications with cross-reactivity to human, mouse, and rat. Suitable for post-transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	HNRNPA0 Rabbit pAb
Remarks/Alias	HNRPA0; HNRNPA0
Species	Human
Gene ID (Human)	10949
Gene ID	10949
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human HNRNPA0 (NP_006796.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13151
Protein Weight	31kDa
Shipping	Ice bag