

Anti-Hnrnpa3 Rabbit Polyclonal Antibody For Wb - P51991

Item Number: CM0021283

Introduction

Rabbit polyclonal antibody targeting human HNRNPA3 (P51991). Validated for WB, IHC-P, ELISA; cross-reacts with mouse, rat. High specificity, ideal for studying RNA trafficking and splicing.

[Learn More](#)

Parameter	Value
Product Name	HNRNPA3 Rabbit pAb
Remarks/Alias	FBRNP; HNRPA3; D10S102; 2610510D13Rik; HNRNPA3
Species	Human
Gene ID	220988
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human HNRNPA3 (NP_001317178.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51991
Protein Weight	40kDa
Shipping	Ice bag