

Anti-Hnrnpcl1 Rabbit Polyclonal Antibody For Wb And Elisa - O60812

Item Number: CM0023030

Introduction

Rabbit polyclonal antibody against human HNRNPCL1, validated for Western blot (WB) and ELISA. Suitable for nucleosome assembly and chromatin research. UniProt O60812, 32kDa.

[Learn More](#)

Parameter	Specification
Product Name	HNRNPCL1 Rabbit pAb
Remarks/Alias	HNRPCL1; HNRNPCL1
Species	Human
Gene ID (Human)	343069
Gene ID	343069
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human HNRNPCL1 (NP_001013653.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O60812
Protein Weight	32kDa
Shipping	Ice bag