

Anti-Hnrnpd Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q14103

Item Number: CM0018142



Introduction

Rabbit polyclonal antibody against HNRNPd (AUF1) for WB, IHC-P, IF/ICC, IP, and ELISA. Validated in human, mouse, and rat samples. Ideal for RNA biology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	HNRNPd Rabbit pAb
Remarks/Alias	P37; AUF1; AUF1A; HNRPD; hnRNPd0; HNRNPd
Species	Human
GeneID (Human)	3184
GeneID	3184
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-306 of human HNRNPd (NP_002129.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14103
Protein Weight	38kDa
Shipping	Ice bag