

Anti-Hnrnpd Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q14103

Item Number: CM0005354



Introduction

Rabbit monoclonal antibody against human HNRNPD (AUF1), validated for WB, IHC-P, IF/ICC, IP, and ELISA. Ideal for mRNA stability and circadian rhythm studies.

[Learn More](#)

Parameter	Value
Product Name	HNRNPD Rabbit mAb
Remarks/Alias	P37; AUF1; AUF1A; HNRPD; hnRNP0; HNRNPD
Species	Human
Gene ID (Human)	3184
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human HNRNPD (NP_112738.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14103
Protein Weight	38 kDa
Shipping	Ice bag