


Anti-Hnrnpd Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14103

Item Number: CM0022908



Introduction

Rabbit polyclonal antibody targeting human HNRNPD (AUF1), validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat samples. Key for studying AU-rich element-mediated mRNA decay and circadian rhythms.

[Learn More](#)

Parameter	Value
Product Name	HNRNPD Rabbit pAb
Remarks/Alias	P37; AUF1; AUF1A; HNRPD; hnRNPD0; HNRNPD
Species	Human
Gene ID	3184 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-306 of human HNRNPD (NP_002129.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14103
Protein Weight	38kDa
Shipping	Ice bag