

Anti-Hnrnpdl Polyclonal Antibody For Wb, If/Icc, Elisa - O14979

Item Number: CM0018560



Introduction

Rabbit polyclonal antibody against human HNRNPDL, validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for transcription and RNA-binding studies.

[Learn More](#)

Parameter	Value
Product Name	HNRNPDL Rabbit pAb
Remarks/Alias	HNRNP; JKTBP; HNRPD; JKTBP2; LGMD1G; LGMD3; IaAUF1; HNRNPDL
Species	Human
Gene ID (Human)	9987
Gene ID	9987
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 261-420 of human HNRNPDL (NP_112740.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14979
Protein Weight	46kDa
Shipping	Ice bag