


Anti-Hnrnp1 Polyclonal Antibody For Wb - P31943

Item Number: CM0030526



Introduction

High-quality rabbit polyclonal antibody targeting human HNRNPH1. Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying pre-mRNA splicing and hnRNP complexes.

[Learn More](#)

Parameter	Value
Product Name	HNRNPH1 Rabbit pAb
Remarks/Alias	HNRPH; HNRPH1; NEDCDS; hnRNP; HNRNPH1
Species	Human
GeneID	3187
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 191-449 of human HNRNPH1 (NP_005511.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31943
Protein Weight	49kDa
Shipping	Ice bag