

Anti-Hnrnpl Polyclonal Antibody For Wb, If/Icc, Elisa - P14866

Item Number: CM0026432



Introduction

Rabbit polyclonal antibody targeting human HNRNPL, an RNA-binding protein involved in splicing and transcriptional regulation. Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for gene expression and neurobiology studies.

[Learn More](#)

Parameter	Value
Product Name	HNRNPL Rabbit pAb
Remarks/Alias	HNRPL; hnRNP-L; P/OKcl.14; HNRNPL
Species	Human
Gene ID (Human)	3191
Gene ID	3191
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 282-589 of human HNRNPL (NP_001524.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14866
Protein Weight	64kDa
Shipping	Ice bag