

## Anti-Hnrrnpm Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P52272

Item Number: CM0030164



### Introduction

High-quality rabbit polyclonal antibody against HNRNPM, validated for WB, IHC-P, IF/ICC, ELISA with cross-reactivity to human, mouse, and rat. Ideal for RNA splicing and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | HNRNPM Rabbit pAb                                                                                                   |
| Remarks/Alias    | CEAR, HNRPM, HTGR1, NAGR1, HNRPM4, HNRNPM4, hnRNP M, HNRNPM                                                         |
| Species          | Human                                                                                                               |
| Gene ID          | 4670                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 491-730 of human HNRNPM (NP_005959.2) |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                            |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | P52272                                                                                                              |
| Protein Weight   | 78kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |