


Anti-Hnrnpr Polyclonal Antibody For Wb, Ihc-P, Elisa - O43390

Item Number: CM0028163



Introduction

High-quality rabbit polyclonal antibody against HNRNPR, validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for mRNA processing and spliceosome studies.

[Learn More](#)

Parameter	Value
Product Name	HNRNPR Rabbit pAb
Remarks/Alias	HNRPR; NEDDFSB; hnRNP-R; HNRNPR
Species	Human
GeneID	10236
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 500 to the C-terminus of human HNRNPR (NP_005817.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43390
Protein Weight	71 kDa
Shipping	Ice bag