


Anti-Hnrnpr Rabbit Polyclonal Antibody [KO Validated] For Wb, If/Icc, Elisa - O43390

Item Number: CM0023448



Introduction

Rabbit polyclonal antibody targeting HNRNPR, a key component of ribonucleosomes involved in pre-mRNA processing. Validated for WB, IF/ICC, ELISA in Human, Mouse, Rat. Ideal for splicing and neurological disease research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HNRNPR Rabbit pAb
Remarks/Alias	HNRPR; NEDDFSB; hnRNP-R; PR
Species	Human
GeneID (human)	10236
GeneID	10236
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 550-633 of human HNRNPR (NP_005817.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43390
Protein Weight	71kDa
Shipping	Ice bag