


Anti-Hnrnpul1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Buj2

Item Number: CM0024231



Introduction

HNRNPUL1 Rabbit Polyclonal Antibody for WB and ELISA. Recognizes human and mouse HNRNPUL1, a transcriptional regulator and RNA-binding protein. Ideal for diagnostic research in transcription and mRNA processing.

[Learn More](#)

Parameter	Value
Product Name	HNRNPUL1 Rabbit pAb
Remarks/Alias	E1BAP5; E1B-AP5; HNRPUL1; HNRNPUL1
Species	Human
GeneID (human)	11100
GeneID	11100
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 209-410 of human HNRNPUL1 (NP_008971.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BUJ2
Protein Weight	42-96 kDa
Shipping	Ice bag