

Anti-Hnrpab Polyclonal Antibody For Wb, If/Icc, Elisa - Q99729

Item Number: CM0021960



Introduction

Rabbit polyclonal antibody against HNRNPAB (HNRNPAB). Validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Recombinant immunogen. For studying RNA binding, hnRNP function, and APOB mRNA editing.

[Learn More](#)

Parameter	Value
Product Name	HNRNPAB Rabbit pAb
Remarks/Alias	ABBP1; HNRNPAB; HNRNPAB
Species	Human
Gene ID (Human)	3182
Gene ID	3182
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-285 of human HNRNPAB (NP_004490.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99729
Protein Weight	36kDa
Shipping	Ice bag