


Anti-Hnrpll Rabbit Polyclonal Antibody For Wb, Elisa - Q8Wvv9

Item Number: CM0018671



Introduction

Anti-HNRPLL Rabbit Polyclonal Antibody validated for WB and ELISA. Recognizes human, mouse, and rat HNRPLL (stromal RNA-regulating factor), a key regulator of alternative splicing in various tissues.

[Learn More](#)

Parameter	Value
Product Name	HNRPLL Rabbit pAb
Remarks/Alias	SRRF; HNRPLL
Species	Human
GeneID (human)	92906
GeneID	92906
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-350 of human HNRPLL (NP_612403.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WVV9
Protein Weight	60kDa
Shipping	Ice bag