

Anti-Homer1 Rabbit Monoclonal Antibody For Wb, Elisa - Q86Ym7

Item Number: CM0003708

Introduction

Recombinant rabbit monoclonal antibody against human Homer1 (HOMER, SYN47). Validated for Western blot (WB) and ELISA, with cross-reactivity to human and rat. Suitable for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	Homer1 Rabbit mAb
Remarks/Alias	HOMER; SYN47; Ves-1; HOMER1A; HOMER1B; HOMER1C; Homer1
Species	Human
GeneID	9456
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 101-220 of human HOMER1 (NP_004263.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q86YM7
Protein Weight	40kDa
Shipping	Ice bag