

Anti-Homer1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q86Ym7

Item Number: CM0010864

Introduction

Rabbit monoclonal antibody against Homer1 (HOMER, SYN47, Ves-1). Validated for WB, IF-P, IHC-P, ELISA in human and rat. Ideal for studying synaptic metabotropic glutamate receptor signaling and neuronal plasticity.

[Learn More](#)

Parameter	Specification
Product Name	Homer1 Rabbit mAb
Remarks/Alias	HOMER; SYN47; Ves-1; HOMER1A; HOMER1B; HOMER1C
Species	Human
GeneID (human)	9456
GeneID	9456
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-354 of human HOMER1 (NP_004263.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q86YM7
Protein Weight	40kDa
Shipping	Ice bag