

Anti-Hooc4 Polyclonal Antibody For Wb, Elisa - P09017

Item Number: CM0028103

Introduction

High-affinity rabbit polyclonal antibody against human HOXC4, validated for Western blot (WB) and ELISA. Reacts with human and mouse HOXC4, a homeobox transcription factor critical for anterior-posterior axis development.

[Learn More](#)

Parameter	Value
Product Name	HOXC4 Rabbit pAb
Remarks/Alias	HOX3, cp19, HOX3E, HOXC4
Species	Human
GeneID (Human)	3221
GeneID	3221
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-130 of human HOXC4 (NP_055435.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P09017
Protein Weight	30kDa
Shipping	Ice bag