

Anti-Hook1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ujc3

Item Number: CM0029596

Introduction

HOOK1 Rabbit Polyclonal Antibody validated for WB and ELISA applications. Reacts with human, mouse, and rat samples. Targets Hook homolog 1, a protein involved in vesicle trafficking and spermatogenesis. Ideal for developmental and cell biology studies.

[Learn More](#)

Parameter	Value
Product Name	HOOK1 Rabbit pAb
Remarks/Alias	HK1; HOOK1
Species	Human
GeneID (Human)	51361
GeneID	51361
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 489-728 of human HOOK1 (NP_056972.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UJC3
Protein Weight	85kDa
Shipping	Ice bag