

Anti-Hook2 Rabbit Monoclonal Antibody For Wb, Elisa - Q96Ed9

Item Number: CM0005970

Introduction

HOOK2 Rabbit Monoclonal Antibody for WB and ELISA applications. Detects human, mouse, and rat HOOK2, a centrosome and vesicle trafficking protein. Suitable for studies of aggresome formation and AP-4 complex distribution.

[Learn More](#)

Parameter	Value
Product Name	HOOK2 Rabbit mAb
Remarks/Alias	HK2; HOOK2
Species	Human
GeneID (Human)	29911
GeneID	29911
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HOOK2 (Q96ED9).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96ED9
Protein Weight	83kDa
Shipping	Ice bag