


Anti-Hook3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q86Vs8

Item Number: CM0013288



Introduction

Rabbit anti-HOOK3 polyclonal antibody validated for WB, IF/ICC, and ELISA. Detects human, mouse, and rat HOOK3, a dynein adapter regulating intracellular trafficking and Golgi organization. Ideal for cytoskeleton and vesicle transport studies.

[Learn More](#)

Parameter	Value
Product Name	HOOK3 Rabbit pAb
Remarks/Alias	HK3; HOOK3
Species	Human
GeneID (Human)	84376
GeneID	84376
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 489-718 of human HOOK3 (NP_115786.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86VS8
Protein Weight	83kDa
Shipping	Ice bag