

Anti-Hoxa11 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P31270

Item Number: CM0030623

Introduction

Anti-HOXA11 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Recognizes human, mouse, rat HOXA11, a homeobox transcription factor essential for anterior-posterior developmental patterning.

[Learn More](#)

Parameter	Value
Product Name	HOXA11 Rabbit pAb
Remarks/Alias	HOX1; HOX1I; RUSAT1; HOXA11
Species	Human
Gene ID (Human)	3207
Gene ID	3207
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human HOXA11 (NP_005514.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31270
Protein Weight	34kDa
Shipping	Ice bag