


Anti-Hoxa13 Monoclonal Antibody For Wb, Elisa - P31271

Item Number: CM0004853



Introduction

HOXA13 rabbit monoclonal antibody validated for WB and ELISA, reacts with human, mouse, and rat. Suitable for studying developmental transcription regulation. SWISS: P31271.

[Learn More](#)

Parameter	Value
Product Name	HOXA13 Rabbit mAb
Remarks/Alias	HOX1; HOX1J; HOXA13
Species	Human
GeneID (human)	3209
GeneID	3209
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 250-350 of human HOXA13 (P31271)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31271
Protein Weight	40kDa
Shipping	Ice bag