


Anti-Hoxa13 Polyclonal Antibody For Wb, Elisa - P31271

Item Number: CM0013749



Introduction

High-quality rabbit polyclonal antibody targeting human HOXA13, validated for Western blot and ELISA. Detects a ~40 kDa band, suitable for developmental biology and transcription factor research.

[Learn More](#)

Parameter	Value
Product Name	HOXA13 Rabbit pAb
Remarks/Alias	HOX1; HOX1J; HOXA13
Species	Human
Gene ID (Human)	3209
Gene ID	3209
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 140-320 of human HOXA13 (NP_000513.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P31271
Protein Weight	40kDa
Shipping	Ice bag