


Anti-Hoxa13 Rabbit Polyclonal Antibody For Wb And Elisa - P31271

Item Number: CM0025827



Introduction

CamelBio offers anti-HOXA13 rabbit polyclonal antibody validated for WB and ELISA, detecting human HOXA13 with cross-reactivity to mouse and rat, facilitating research on developmental transcription factors and positional identity.

[Learn More](#)

Parameter	Value
Product Name	HOXA13 Rabbit pAb
Remarks/Alias	HOX1; HOX1J; HOXA13
Species	Human
GeneID (Human)	3209
GeneID	3209
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 140-320 of human HOXA13 (NP_000513.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P31271
Protein Weight	40kDa
Shipping	Ice bag