

Anti-Hoxa7 Rabbit Polyclonal Antibody For Wb, Elisa - P31268

Item Number: CM0017147



Introduction

High-quality rabbit polyclonal antibody targeting human HOXA7, validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for studying developmental gene regulation and anterior-posterior axis formation.

[Learn More](#)

Parameter	Value
Product Name	HOXA7 Rabbit pAb
Remarks/Alias	ANTP; HOX1; HOX1A; HOX1.1; HOXA7
Species	Human
Gene ID (Human)	3204
Gene ID	3204
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human HOXA7 (NP_008827.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31268
Protein Weight	25kDa
Shipping	Ice bag