

Anti-Hoxa9 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P31269

Item Number: CM0002095

Introduction

High-quality rabbit monoclonal antibody against human HOXA9 (homeobox protein Hox-A9, P31269), validated for WB, IHC-P, and ELISA. Cross-reacts with rat. Essential for developmental biology, oncology, and transcription research.

[Learn More](#)

Parameter	Value
Product Name	HOXA9 Rabbit mAb
Remarks/Alias	HOX1; ABD-B; HOX1G; HOX1.7; HOXA9
Species	Human
Gene ID (Human)	3205
Gene ID	3205
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 173-272 of human HOXA9 (P31269).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P31269
Protein Weight	30kDa
Shipping	Ice bag