

Anti-Hoxa9 Rabbit Polyclonal Antibody For Wb And Elisa - P31269

Item Number: CM0022422

Introduction

Rabbit polyclonal antibody targeting human HOXA9 homeobox transcription factor. Validated for WB and ELISA, reacts with human and mouse. Suitable for developmental biology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	HOXA9 Rabbit pAb
Remarks/Alias	HOX1; ABD-B; HOX1G; HOX1.7; HOXA9
Species	Human
Gene ID (Human)	3205
Gene ID	3205
Immunogen	Recombinant protein matching amino acids 1-190 of human HOXA9 (NP_689952.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P31269
Protein Weight	30kDa
Shipping	Ice bag