

Anti-Hoxb1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P14653

Item Number: CM0027547



Introduction

Rabbit polyclonal antibody against human HOXB1, validated for Western blot, immunofluorescence/immunocytochemistry, and ELISA. Reacts with human, mouse, and rat. Suitable for developmental biology and transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	HOXB1 Rabbit pAb
Remarks/Alias	HOX2, HCFP3, HOX2L, Hox-2.9, HOXB1
Species	Human
Gene ID (Human)	3211
Gene ID	3211
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14653
Protein Weight	32kDa
Shipping	Ice bag