

Anti-Hoxb2 Rabbit Polyclonal Antibody For Wb And Elisa - P14652

Item Number: CM0021556

Introduction

Anti-HOXB2 Rabbit Polyclonal Antibody validated for Western blot (WB) and ELISA applications. Detects human Homeobox protein Hox-B2 and cross-reacts with mouse and rat. Ideal for studies of anterior-posterior axis development.

[Learn More](#)

Parameter	Value
Product Name	HOXB2 Rabbit Polyclonal Antibody
Remarks/Alias	K8; HOX2; HOX2H; Hox-2.8; HOXB2
Species	Human
Gene ID	3212
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-140 of human HOXB2 (NP_002136.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P14652
Protein Weight	38kDa
Shipping	Ice bag