

Anti-Hoxb4 Rabbit Polyclonal Antibody For Wb, Elisa - P17483

Item Number: CM0024240

Introduction

Anti-HOXB4 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human and mouse HOXB4. Ideal for studying developmental transcription regulation and anterior-posterior axis patterning.

[Learn More](#)

Parameter	Value
Product Name	HOXB4 Rabbit pAb
Remarks/Alias	HOX2; HOX2F; HOX-2.6; HOXB4
Species	Human
Gene ID (Human)	3214
Gene ID	3214
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human HOXB4 (NP_076920.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P17483
Protein Weight	28kDa
Shipping	Ice bag