

Anti-Hoxb6 Polyclonal Antibody For Wb, Elisa - P17509

Item Number: CM0019380

Introduction

Rabbit polyclonal antibody targeting human HOXB6, a homeobox transcription factor. Validated for Western blot and ELISA, cross-reacts with mouse. Ideal for developmental and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	HOXB6 Rabbit pAb
Remarks/Alias	HOX2; HU-2; HOX2B; Hox-2.2; HOXB6
Species	Human
Gene ID	3216
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-110 of human HOXB6 (NP_061825.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P17509
Protein Weight	25kDa
Shipping	Ice bag