

## Anti-Hoxb7 Rabbit Polyclonal Antibody For Wb, Elisa - P09629

Item Number: CM0027635

### Introduction

Rabbit polyclonal antibody against human HOXB7, validated for WB and ELISA. Detects human and mouse HOXB7. Suitable for developmental biology and cancer research. Immunogen: N-terminal recombinant protein.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | HOXB7 Rabbit pAb                                                                                                  |
| Remarks/Alias    | HOX2; HOX2C; HHO.C1; Hox-2.3; HOXB7                                                                               |
| Species          | Human                                                                                                             |
| Gene ID (Human)  | 3217                                                                                                              |
| Gene ID          | 3217                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human HOXB7 (NP_004493.3). |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse                                                                                                      |
| SWISS            | P09629                                                                                                            |
| Protein Weight   | 24kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |